
NESS L300[®] Plus

NESS L300[®] Plus System

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Manchon de cuisse

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NESS L300[®]Plus



Intelli-Gait Plus Software

© 2010 Bioness Neuromodulation Ltd., a Bioness Inc. Company. All rights reserved

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Phase d'appui

The screenshot displays the 'Gait Settings' interface for a patient named John D. The interface includes a top status bar with a Windows logo, the title 'Gait Settings - John D.', and a timer showing '03:00'. Below this, there are two muscle selection options: 'Thigh' and 'Hams.'. The 'Hams.' option is currently selected, with a value of '47' displayed. A 'Start' button is visible. The interface also features sliders for 'Delay (%)' and 'End (%)', both set to 10. A central graph shows a gait cycle with 'Swing Control' and 'Stance Control' checkboxes. The 'Stance Control' checkbox is checked. A 'Default' button is also present. At the bottom, there is a navigation bar with 'Exit', 'Patients', 'Settings', 'History', and 'Tools'.

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Default

Setup Stim Gait Training

Exit | Patients | Settings | History | Tools

Paramètres de marche avancés, ischios-jambiers

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The screenshot displays the 'Adv. Gait Settings' interface for a patient named 'John D.'. The interface is organized into several sections:

- Header:** A blue bar at the top contains the Windows logo, the text 'Adv. Gait Settings - John D.', a refresh icon, a volume icon, and the time '11:33'.
- Muscle Selection:** On the left, there are two anatomical diagrams of a leg. The top one is labeled 'Thigh' and shows the quadriceps muscle. The bottom one is labeled 'L300' and shows the hamstring muscle. A vertical slider between them is labeled 'Hams.' and has the number '51' in the center, indicating the current selection level.
- Stimulation Parameters:** In the center, there are four trapezoidal graphs representing different phases of the gait cycle: 'Swing Ramp Up', 'Swing Ramp Down', 'Stance Ramp Up', and 'Stance Ramp Down'. Each graph has a numerical value '0.1' in a box with left and right arrow buttons, indicating the ramp rate.
- Max. Duration of Stimulation:** Below the graphs, there is a section titled 'Max. Duration of Stimulation' with a numerical input field showing '4' and the unit 'Seconds'. A 'Default' button is located to the right of this section.
- Navigation:** At the bottom, there are four large buttons labeled 'Setup', 'Stim', 'Gait', and 'Training'. Below these is a light blue bar with navigation links: 'Exit | Patients | Settings | History | Tools'.
- Start Button:** A prominent green 'Start' button with a green indicator light is located in the upper right area of the main control panel.

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Ramp Down 1 (0 - 2 Sec)

Total Time 20 (5-60 Min)

Enable specific training intensity ?

Start

Both

Thigh Only

L300 Only

Default

Setup Stim Gait Training

Exit | Patients | Settings | History | Tools

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
 **Gait Log - John D.**   **01:28**

From: To: Resolution:  


Date	# of Steps	Total	Gait Duration	Total
06/11/11	0	0	00:00:00	00:00:00
06/12/11	0	0	00:00:00	00:00:00
06/13/11	58	58	00:00:54	00:00:54
06/14/11	144	202	00:02:35	00:03:29
06/15/11	276	478	00:03:57	00:07:26
06/16/11	407	885	00:04:39	00:12:05
06/17/11	556	1441	00:08:52	00:20:57
Average:	206		00:03:00	



Exit | **Patients** | **Settings** | **History** | **Tools** 

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